

Flask. (Save a little liquid in pipette each time.)

Tips: A 10 mL pipette allows media to be dispensed at the bottom of the vessel. A 25 mL pipette allows media to be dispensed just past the NEST Logo.

3. Hold the Multi-layer Flask upright with the Logo facing you and tilt clockwise to a 45° angle on a flat work surface to partition the liquid evenly into each layer.

4. While holding the Multi-layer Flask at a 45° angle, gently lay it flat onto the work surface with logo facing up.

5. After placing the Multi-layer Flask flat on a work surface, gently rock back and forth and side-to-side to distribute cells evenly onto culture surfaces.

6. Put the flask quickly and slightly into the incubator.

Tips: Be careful to avoid foaming of medium, and not to spill liquid from each layer.

Removal of culture medium:

You may choose to either aspirate or pour the media from Multi-layer Flask.

Aspirating method: To aspirate or remove media, Hold the Multi-layer Flask upright, with the NEST Logo facing you, tilt counter-clock wise to a 45° angle. Then, Extend the pipette to the bottom to aspirate the medium until it is completely aspirated.

Pouring method: Hold the Multi-layer Flask upright, with Logo facing you, tilt counter-clock wise to a 45° angle and pour spent media from Multi-layer Flask

Tips: It is recommended to use a NEST 10mL aspirating pipette when aspirating the culture medium, so that the medium can be completely sucked up.

Collection of cells:

1. Wash with buffer for one time and add dissociating reagent (≥5mL per layer). Then, follow Steps 3-5 to distribute to dissociating reagent to each layer.

2. Neutralize with inactivating solution and mix following Steps 3-5. Gently swirl to dislodge cells completely.

3. Follow "Aspirating Method" protocol and collect cell suspension using a NEST 10mL serological pipette.

4. Follow "Pouring Method". Pour the cell suspension into a NEST conical tube.

5. Rinse with additional wash buffer if needed.

6. Search "NEST Multi-layer Flask" video on NEST website or Youtube.

Recommendations

Please be gently when use the 5-layer cell culture flask to avoid air bubbles. Air bubbles will lead the culture medium flow from the upper layer to the lower layer.

Product Range

3-layer Cell Culture Flask


Cat. No.	Growth area	Cap Style	Size (mm)			TC Treated	/Pack	/Case
			Length	Width	Height			
CS0007-731301	520 cm ²	Plug Seal Cap	203	120.5	60.1	Yes	1	12
CS0007-731302		Vent Cap						

5-layer Cell Culture Flask

Cat. No.	Growth area	Cap Style	Size (mm)			TC Treated	/Pack	/Case
			Length	Width	Height			
CS0007-731001	870 cm ²	Plug Seal Cap	203.6	120.5	84.3	Yes	1	8
CS0007-731002		Vent Cap						

Technical Drawing

3-layer Cell Culture Flask



5-layer Cell Culture Flask

