

888-466-3604 | support@bosterbio.com | www.bosterbio.com

# Anti-iNOS/NOS2 Antibody Picoband<sup>™</sup>

Catalog Number: A00368

## About NOS2

Nitric oxide synthase, inducible is an enzyme that in humans is encoded by the NOS2 gene. Nitric oxide (NO) is a messenger molecule with diverse functions throughout the body. In the brain and peripheral nervous system, NO displays many properties of a neurotransmitter; it is implicated in neurotoxicity associated with stroke and neurodegenerative diseases, neural regulation of smooth muscle, including peristalsis, and penile erection. Three different NOS isoforms have been identified which fall into two distinct types, constitutive and inducible. The inducible NOS (iNOS) isoform is expressed in a variety of cell types and tissues in response to inflammatory agents and cytokines. The human iNOS (NOS2) gene is approximately 37 kb in length and consists of 26 exons and 25 introns. NOS2-derived NO is a prerequisite for cytokine signaling and function in innate immunity.

### **Overview**

Product Name	Anti-iNOS/NOS2 Antibody Picoband™
Reactive Species	Human
Description	Boster Bio Anti-iNOS/NOS2 Antibody Picoband <sup>™</sup> catalog # A00368. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains 5mg BSA, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35228

### **Technical Details**

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human iNOS, different from the related mouse sequence by five amino acids, and from the related rat sequence by four amino acids.
Predicted Reactive Species	Chicken
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG



# BOSTER BIOLOGICAL TECHNOLOGY 3942 B Valley Ave, Pleasanton, CA 94566

888-466-3604 | support@bosterbio.com | www.bosterbio.com

Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: Western blot, 0.1-0.5ug/ml, Human



888-466-3604 | support@bosterbio.com | www.bosterbio.com

### Anti-iNOS/NOS2 Antibody Picoband<sup>™</sup> (A00368) Images

250KD -	
	Figure 1. Western blot analysis of iNOS using anti-iNOS
130KD	antibody (A00368). Electrophoresis was performed on a 5-20% SDS-PAGE gel at
100KD -	70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The
	sample well of each lane was loaded with 50ug of sample
70KD -	under reducing conditions.
55KD -	Lane 1: HELA whole cell lysates.
SOLD	After Electrophoresis, proteins were transferred to a
	Nitrocellulose membrane at 150mA for 50-90 minutes.
	Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5
	hour at RT. The membrane was incubated with rabbit anti-
	iNOS antigen affinity purified polyclonal antibody (Catalog #
	A00368) at 0.5 ug/mL overnight at 4°C, then washed with
	TBS-0.1%Tween 3 times with 5 minutes each and probed
	with a goat anti-rabbit IgG-HRP secondary antibody at a
	dilution of 1:10000 for 1.5 hour at RT. The signal is
	developed using an Enhanced Chemiluminescent detection
	(ECL) kit (Catalog # EK1002) with Tanon 5200 system. A
	specific band was detected for iNOS at approximately
	130KD. The expected band size for iNOS is at 130KD.

### **34 Publications Citing This Product**

1. PubMed ID: 10.3389/fphar.2017.00044, Celastrol Attenuates Multiple Sclerosis and Optic Neuritis in an Experimental Autoimmune Encephalomyelitis Model

2. PubMed ID: 10.3109/13880209.2014.898075, Effect of ammonium pyrrolidine dithiocarbamate (PDTC) on NF-kappaB activation and CYP2E1 content of rats with immunological liver injury

3. PubMed ID: 10.3892/mmr.2014.2859, HIF21 signaling pathway involving iNOS, COX22 and caspase29 mediates the neuroprotection provided by erythropoietin in the retina of chronic ocular hypertension rats

Visit bosterbio.com/anti-inos-picoband-trade-antibody-a00368-boster.html to see all 34 publications.

### Submit a product review to Biocompare.com



Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-iNOS/NOS2 Antibody ™