

Anti-Angiotensinogen AGT Antibody

Catalog Number: A02103-2

About AGT

Chromatin reader that specifically recognizes and binds histone H3.3 trimethylated at 'Lys-36' (H3.3K36me3) and regulates RNA polymerase II elongation. Does not bind other histone H3 subtypes (H3.1 or H3.2) By similarity. Colocalizes with highly expressed genes and functions as a transcription corepressor by modulating RNA polymerase II at the elongation stage. Acts as a tumor-suppressor by repressing a transcriptional program essential for tumor cell growth.

Hateboer G., EMBO J. 14:3159-3169(1995).

Deloukas P., Nature 429:375-381(2004).

Matsuoka S., Science 316:1160-1166(2007).

Overview

Product Name	Anti-Angiotensinogen AGT Antibody
Reactive Species	Human
Description	Boster Bio Anti-Angiotensinogen AGT Antibody catalog # A02103-2. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01019

Technical Details

Immunogen	Synthesized peptide derived from human Siah-1/2
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified

Suggested Dilutions

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

WB: 1:1,000

ICC: 1:50-1:200

Anti-Angiotensinogen AGT Antibody (A02103-2) Images




Figure 1. Immunocytochemistry staining of AGT using Anti-Angiotensinogen AGT Antibody (A02103-2).

ICC staining Angiotensinogen in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 2. Immunocytochemistry staining of AGT using Anti-Angiotensinogen AGT Antibody (A02103-2).

ICC staining Angiotensinogen in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Angiotensinogen AGT Antibody