

Anti-IA-2/PTPRN Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A03116-1-Fluoro647

About PTPRN

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and a single catalytic domain, and thus represents a receptor-type PTP. This PTP was found to be an autoantigen that is reactive with insulin-dependent diabetes mellitus (IDDM) patient sera, and thus may be a potential target of autoimmunity in diabetes mellitus. Alternate splicing results in multiple transcript variants.

Overview

| | |
|----------------------|--|
| Product Name | Anti-IA-2/PTPRN Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q16849 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human IA-2/PTPRN recombinant protein (Position: V35-H569). Human IA-2/PTPRN shares 79.5% amino acid (aa) sequence identity with rat IA-2/PTPRN. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IA-2/PTPRN Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.