

Anti-PDE8A Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05304-1-Fluoro488

About PDE8A

High affinity cAMP-specific and IBMX-insensitive 3', 5'-cyclic phosphodiesterase 8A is an enzyme that in humans is encoded by the PDE8A gene. The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase (PDE) family, and PDE8 subfamily. This PDE hydrolyzes the second messenger, cAMP, which is a regulator and mediator of a number of cellular responses to extracellular signals. Thus, by regulating the cellular concentration of cAMP, this protein plays a key role in many important physiological processes. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-PDE8A Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O60658

Technical Details

Immunogen	E.coli-derived human PDE8A recombinant protein (Position: R32-E495). Human PDE8A shares 80.5% amino acid (aa) sequence identity with mouse PDE8A.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDE8A Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.