

Anti-NFKBIZ Antibody Picoband® Cy3 Conjugated

Catalog Number: A05312-1-Cy3

About NFKBIZ

NF-kappa-B inhibitor zeta (IkappaBzeta) is a protein that in humans is encoded by the NFKBIZ gene. This gene is a member of the ankyrin-repeat family and is induced by lipopolysaccharide (LPS). The C-terminal portion of the encoded product which contains the ankyrin repeats, shares high sequence similarity with the I kappa B family of proteins. The latter are known to play a role in inflammatory responses to LPS by their interaction with NF-B proteins through ankyrin-repeat domains. Studies in mouse indicate that this gene product is one of the nuclear I kappa B proteins and an activator of IL-6 production. Two transcript variants encoding different isoforms have been found for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-NFKBIZ Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9BYH8 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human NFKBIZ recombinant protein (Position: Q247-H577). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NFKBIZ Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.