



Wash Buffer (25x)

Catalog number: AR1106-5

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Wash Buffer (25x)

Catalog Number: AR1106-5, **Storage:** Upon receipt store at 2-8°C. It is stable at -2-8°C for 24 months.

Product Name	Wash Buffer (25x)
SKU/Catalog Number	AR1106-5
Pack Size	20 ml
Storage	Upon receipt store at 2-8°C. It is stable at -2-8°C for 24 months.
Usage Method	Prepare of Working Wash Buffer by diluting the supplied 20 ml of Wash Buffer (25x) with 480ml of deionized or water. If crystals have formed in the concentrate, warm to room temperature and mix it gently until crystals have completely dissolved.
Cite This Product	Wash Buffer (25x) (Boster Biological Technology, Pleasanton CA, USA, Catalog # AR1106-5)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Wash Buffer (25x)

For Research Use Only. Not for use in diagnostic procedures.