

Anti-TNFSF2/TNFA Reference Antibody (adalimumab)

Catalog Number: M00002-5


Overview

| | |
|----------------------|--|
| Product Name | Anti-TNFSF2/TNFA Reference Antibody (adalimumab) |
| Reactive Species | Cynomolgus monkey, Human, Mouse |
| Description | Boster Bio Anti-TNFSF2/TNFA Reference Antibody (adalimumab) (Catalog # M00002-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.001EU/ug,determined by LAL method. Expression system: CHO Cell |
| Conjugate | Unconjugated |
| Application | ELISA, Flow Cytometry, Functional Assay, Kinetics |
| Clonality | Monoclonal |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles. |
| Uniprot ID | P01375 |


Technical Details

| | |
|------|---|
| Form | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
|------|---|


Anti-TNFSF2/TNF α Reference Antibody (adalimumab) (M00002-5) Images




Anti-TNFSF2 / TNF α Reference Antibody (adalimumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



Anti-TNFSF2/TNF α Reference Antibody (adalimumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFSF2/TNF α Reference Antibody (adalimumab) is more than 98.86%



The purity of Anti-TNFSF2 / TNF α Reference Antibody (adalimumab) is more than 98.86% ,determined by SEC-HPLC.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFSF2/TNF α Reference Antibody (adalimumab)

For Research Use Only. Not for use in diagnostic procedures.