

Anti-p53 (phospho-S20) TP53 Antibody

Catalog Number: A00001S20


Overview

Product Name	Anti-p53 (phospho-S20) TP53 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-p53 (phospho-S20) TP53 Antibody catalog # A00001S20. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04637

Technical Details

Immunogen	Synthetic phosphopeptide derived from human p53 around the phosphorylation site of Serine 20.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-p53 (phospho-S20) TP53 Antibody (A00001S20) Images



Immunohistochemistry (IHC) analyzes of p-P53 (S20) pAb in paraffin-embedded human breast carcinoma tissue.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p53 (phospho-S20) TP53 Antibody

For Research Use Only. Not for use in diagnostic procedures.