

## Anti-PTEN Antibody Picoband® Cy3 Conjugated

Catalog Number: A00006-1-Cy3

### About PTEN

PTEN is also known as BZS, DEC, CWS1, GLM2, MHAM, TEP1, PTEN1. It is mapped to 10q23.3. This gene was identified as a tumor suppressor that is mutated in a large number of cancers at high frequency. The protein encoded by this gene is a phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase. It contains a tensin like domain as well as a catalytic domain similar to that of the dual specificity protein tyrosine phosphatases. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. The PTEN structure reveals a phosphatase domain that is similar to protein phosphatases but also has an enlarged active site important for the accommodation of the phosphoinositide substrate.

### Overview

Product Name	Anti-PTEN Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P60484

### Technical Details

Immunogen	E. coli-derived human PTEN recombinant protein (Position: T2-T202). Human PTEN shares 100% amino acid (aa) sequence identity with mouse PTEN.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 14 Publications Citing This Product

1. PubMed ID: 10.3892/ijo.2016.3401, miR-223 reverses the resistance of EGFR-TKIs through IGF1R/PI3K/Akt signaling pathway
2. PubMed ID: 31614022, Sun M,Hu L, S,Huang T,Zhang M,Yang M,Zhen W,Yang D,Lu W,Guan M,Peng S. Circulating MicroRNA-19b Identified From Osteoporotic Vertebral Compression Fracture Patients Increases Bone Formation. J Bone Miner Res.2020 Feb;35(2):306-316. doi:10.1002/jbmr.3892.
3. PubMed ID: 26936292, miR-223 reverses the resistance of EGFR-TKIs through IGF1R/PI3K/Akt signaling pathway

Visit [bosterbio.com/anti-pten-antibody-a00006-1-boster.html](http://bosterbio.com/anti-pten-antibody-a00006-1-boster.html) to see all 14 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PTEN Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.