

## Anti-PTEN Antibody Picoband® HRP Conjugated

Catalog Number: A00006-1-HRP

### About PTEN

PTEN is also known as BZS, DEC, CWS1, GLM2, MHAM, TEP1, PTEN1. It is mapped to 10q23.3. This gene was identified as a tumor suppressor that is mutated in a large number of cancers at high frequency. The protein encoded by this gene is a phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase. It contains a tensin like domain as well as a catalytic domain similar to that of the dual specificity protein tyrosine phosphatases. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. The PTEN structure reveals a phosphatase domain that is similar to protein phosphatases but also has an enlarged active site important for the accommodation of the phosphoinositide substrate.

### Overview

Product Name	Anti-PTEN Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P60484

### Technical Details

Immunogen	E. coli-derived human PTEN recombinant protein (Position: T2-T202). Human PTEN shares 100% amino acid (aa) sequence identity with mouse PTEN.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## 14 Publications Citing This Product

---

1. PubMed ID: 10.3892/ijo.2016.3401, miR-223 reverses the resistance of EGFR-TKIs through IGF1R/PI3K/Akt signaling pathway
2. PubMed ID: 31614022, Sun M,Hu L, S,Huang T,Zhang M,Yang M,Zhen W,Yang D,Lu W,Guan M,Peng S. Circulating MicroRNA-19b Identified From Osteoporotic Vertebral Compression Fracture Patients Increases Bone Formation. J Bone Miner Res.2020 Feb;35(2):306-316. doi:10.1002/jbmr.3892.
3. PubMed ID: 26936292, miR-223 reverses the resistance of EGFR-TKIs through IGF1R/PI3K/Akt signaling pathway

Visit [bosterbio.com/anti-pten-antibody-a00006-1-boster.html](https://bosterbio.com/anti-pten-antibody-a00006-1-boster.html) to see all 14 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PTEN Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.