

Anti-Stat3 (phospho-Y705) Antibody

Catalog Number: A00007Y705-1


Overview

Product Name	Anti-Stat3 (phospho-Y705) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Stat3 (phospho-Y705) Antibody catalog # A00007Y705-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P40763


Technical Details

Immunogen	Synthetic phosphopeptide derived from human Stat3 around the phosphorylation site of Tyrosine 705.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

Anti-Stat3 (phospho-Y705) Antibody (A00007Y705-1) Images



Western blot (WB) analysis of Stat3 (phospho-Y705) polyclonal antibody at 1:500 dilution Lane1:HeLa whole cell lysate(40ug) Lane2:HeLa treated with IFN-gamma(100ng/ml, 30min) whole cell lysate(40ug) Lane3:Jurkat whole cell lysate(40ug) Lane4:Jurkat treated with H2O2(100nM/ml, 30min) whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-STAT3 (Y705) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat3 (phospho-Y705) Antibody

For Research Use Only. Not for use in diagnostic procedures.