

## Anti-STAT3 (Phospho-Tyr705) Antibody

Catalog Number: A00007Y705-2

### Overview


Product Name	Anti-STAT3 (Phospho-Tyr705) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-STAT3 (Phospho-Tyr705) Antibody catalog # A00007Y705-2. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P40763

### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Tyr705 of human STAT3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

## Anti-STAT3 (Phospho-Tyr705) Antibody (A00007Y705-2) Images

---



Immunohistochemical analysis of paraffin-embedded human lung tissue using anti-Phospho-STAT3(Tyr705) antibody. Counter stained with hematoxylin.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STAT3 (Phospho-Tyr705) Antibody

For Research Use Only. Not for use in diagnostic procedures.