

## Anti-Forkhead box protein P3 FOXP3 Antibody

Catalog Number: A00011-1


### Overview

Product Name	Anti-Forkhead box protein P3 FOXP3 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Forkhead box protein P3 FOXP3 Antibody catalog # A00011-1. Tested in ICC, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BZS1


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	ICC: 1:50-1:200 FC: 1:50-1:100

## Anti-Forkhead box protein P3 FOXP3 Antibody (A00011-1) Images



ICC staining FOXP3 in MCF-7 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining FOXP3 in 293T cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

## 5 Publications Citing This Product

1. PubMed ID: 27698748, IL-35 improves Treg-mediated immune suppression in atherosclerotic mice
2. PubMed ID: 26064167, Gong T, Wang Cf, Yuan Jr, Li Y, Gu Jf, Zhao Bj, Zhang L, Jia Xb, Feng L, Liu Sl. Evid Based Complement Alternat Med. 2015;2015:630760. Doi: 10.1155/2015/630760. Epub 2015 May 3. Inhibition Of Tumor Growth And Immunomodulatory Effects Of Flavonoids...
3. PubMed ID: 25197340, Hu M, Li K, Maskey N, Xu Z, Peng C, Wang B, Li Y, Yang G. Int J Clin Exp Pathol. 2014 Jul 15;7(8):4685-94. Ecollection 2014. Decreased Intratumoral Foxp3 Tregs And Increased Dendritic Cell Density By Neoadjuvant Chemotherapy Associated With Favora...

Visit [bosterbio.com/anti-forkhead-box-protein-p3-foxp3-antibody-a00011-1-boster.html](http://bosterbio.com/anti-forkhead-box-protein-p3-foxp3-antibody-a00011-1-boster.html) to see all 5 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Forkhead box protein P3 FOXP3 Antibody

For Research Use Only. Not for use in diagnostic procedures.