

## Anti-HIF-1 alpha/HIF1A Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A00013-3-Fluoro488

### About HIF1A

HIF-1alpha (Hypoxia-inducible factor 1alpha, HIF1A) is a transcription factor that mediates cellular and systemic homeostatic responses to reduced O<sub>2</sub> availability in mammals, including angiogenesis, erythropoiesis and glycolysis. This gene was mapped to 14q21-q24. HIF-1alpha transactivate genes required for energy metabolism and tissue perfusion and is necessary for embryonic development and tumor explant growth. HIF-1alpha is over expressed during carcinogenesis, myocardial infarction and wound healing. It is crucial for the cellular response to hypoxia and is frequently over expressed in human cancers, resulting in the activation of genes essential for cell survival. HIF-1alpha regulates the survival and function in the inflammatory microenvironment directly. It is a transcription factor that plays a pivotal role in cellular adaptation to changes in oxygen availability.

### Overview

Product Name	Anti-HIF-1 alpha/HIF1A Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q16665

### Technical Details

Immunogen	E.coli-derived human HIF-1 alpha/HIF1A recombinant protein (Position: Q603-Q814).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 28 Publications Citing This Product

1. PubMed ID: 10.1590/1414-431X20142910, Putative role of ischemic postconditioning in a rat model of limb ischemia and reperfusion: involvement of hypoxia-inducible factor-1alpha expression
2. PubMed ID: 10.3748/wjg.v21.i7.2191, Semaphorin 4D and hypoxia-inducible factor-1alpha overexpression is related to prognosis in colorectal carcinoma
3. PubMed ID: 10.1186/s13046-016-0399-7, Hyperthermia induced HIF-1a expression of lung cancer through AKT and ERK signaling pathways

Visit [bosterbio.com/anti-hif-1-alpha-hif1a-picoband-trade-antibody-a00013-3-boster.html](http://bosterbio.com/anti-hif-1-alpha-hif1a-picoband-trade-antibody-a00013-3-boster.html) to see all 28 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HIF-1 alpha/HIF1A Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.