

## Anti-TGF Beta 1 (Y369) TGFB1 Antibody

Catalog Number: A00019Y369


### Overview

Product Name	Anti-TGF Beta 1 (Y369) TGFB1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TGF Beta 1 (Y369) TGFB1 Antibody catalog # A00019Y369. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01137

### Technical Details

Immunogen	Synthetic peptide, corresponding to Human TGFbeta1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:5000-1:10000

## Anti-TGF Beta 1 (Y369) TGFB1 Antibody (A00019Y369) Images



Western blot (WB) analysis of TGFbeta1 (Y369) pAb at 1:5000 dilution Lane1:C6 whole cell lysate(40ug) Lane2:BV2 whole cell lysate(40ug) Lane3:A549 whole cell lysate(40ug) Lane4:Hela whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TGF Beta 1 (Y369) TGFB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.