

Anti-EGFR (Phospho-Y1173) Antibody

Catalog Number: A00023Y1173


Overview

Product Name	Anti-EGFR (Phospho-Y1173) Antibody
Reactive Species	Human
Description	Boster Bio Anti-EGFR (Phospho-Y1173) Antibody catalog # A00023Y1173. Tested in WB,ICC,IHC,IP applications. This antibody reacts with Human.
Application	IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P00533


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Tyr1173 of human EGFR.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-EGFR (Phospho-Y1173) Antibody (A00023Y1173) Images



ICC staining (pY1173) in untreated A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining Phospho-EGFR (pY1173) in untreated A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EGFR (Phospho-Y1173) Antibody

For Research Use Only. Not for use in diagnostic procedures.