

Anti-c-Myc (Phospho-T58) Antibody

Catalog Number: A00026T58-1


Overview

Product Name	Anti-c-Myc (Phospho-T58) Antibody
Reactive Species	Human
Description	Boster Bio Anti-c-Myc (Phospho-T58) Antibody catalog # A00026T58-1. Tested in WB,ICC/IF,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01106


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr58 of human c-Myc
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 FC: 1:50-1:100 ICC: 1:50-1:200

Anti-c-Myc (Phospho-T58) Antibody (A00026T58-1) Images



ICC staining phospho-c-Myc(T58) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining phospho-c-Myc(T58) in SKOV-3 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-c-Myc (Phospho-T58) Antibody

For Research Use Only. Not for use in diagnostic procedures.