

## Anti-c-Myc (Phospho-T58/S62) Antibody

Catalog Number: A00026T58S62


### Overview

Product Name	Anti-c-Myc (Phospho-T58/S62) Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-c-Myc (Phospho-T58/S62) Antibody catalog # A00026T58S62. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human,Rat.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01106

### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr58 and Ser62 of human c-Myc.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:50-1:200

## Anti-c-Myc (Phospho-T58/S62) Antibody (A00026T58S62) Images



Western blot analysis of Phospho-c-Myc(T58+S62) on K562 cells lysates using anti-Phospho-c-Myc(T58+S62) antibody at 1/1,000 dilution.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-c-Myc (Phospho-T58/S62) Antibody

For Research Use Only. Not for use in diagnostic procedures.