

Anti-Stat1 (phospho-Y701) Antibody

Catalog Number: A00036Y701-2


Overview

Product Name	Anti-Stat1 (phospho-Y701) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Stat1 (phospho-Y701) Antibody catalog # A00036Y701-2. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42224


Technical Details

Immunogen	Synthetic phosphopeptide derived from human Stat1 around the phosphorylation site of Tyrosine 701.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

Anti-Stat1 (phospho-Y701) Antibody (A00036Y701-2) Images



Western blot (WB) analysis of p-Stat1 (Y701) pAb at 1:500 dilution Lane1:U-87MG whole cell lysate(40ug) Lane2:Hela whole cell lysate(40ug) Lane3:PC12 whole cell lysate(40ug) Lane4:The Embryo tissue lysate of Mouse(40ug) Lane5:Panc1 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-Stat1 (Y701) pAb in paraffin-embedded human breast carcinoma tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat1 (phospho-Y701) Antibody

For Research Use Only. Not for use in diagnostic procedures.