

## Anti-AAG/MPG Antibody Picoband®

Catalog Number: A00037

### About MPG

MPG(N-methylpurine-DNA glycosylase) also known as MDG, 3-METHYLADENINE DNA GLYCOSYLASE, 3MeAde DNA GLYCOSYLASE, AAG or APNG. The MPG gene is mapped to human chromosome 16 by analysis of a panel of DNAs from mouse/human and hamster/human hybrid cell lines. The MPG gene was expressed in all cell lines and tissues examined, but was found at particularly high levels in a colon adenocarcinoma cell line (HT29). The completely characterized human MPG gene was found to span 7 to 8 kb of genomic DNA and to be localized 75 kb upstream of the embryonic zeta-globin gene. To assess the contribution of Apng to the repair of several mutagenic lesions in vivo, Hang et al. (1997) biochemically analyzed cell-free extracts of tissues from mice with a targeted deletion of the Apng gene. Following treatment with DNA-methylating agents, increased persistence of 7-methylguanine (meG) was found in liver sections of APNG knockout mice in comparison with wildtype mice, demonstrating an in vivo phenotype for the APNG-null animals.

### Overview


Product Name	Anti-AAG/MPG Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-AAG/MPG Antibody Picoband® catalog # A00037. Tested in ELISA, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P29372

### Technical Details

Immunogen	E.coli-derived human AAG/MPG recombinant protein (Position: Q20-A298).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG

Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml/ml, Human ELISA, 0.1-0.5 ug/ml/ml, Human

## Anti-AAG/MPG Antibody Picoband® (A00037) Images



Western blot analysis of AAG/MPG using anti-AAG/MPG antibody (A00037). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human HepG2 whole cell lysates, Lane 2: human Jurkat whole cell lysates, Lane 3: human 293T whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-AAG/MPG antigen affinity purified polyclonal antibody (Catalog # A00037) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for AAG/MPG at approximately 37 kDa. The expected band size for AAG/MPG is at 33 kDa.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AAG/MPG Antibody

For Research Use Only. Not for use in diagnostic procedures.