

Anti-Caspase-8 (S347) CASP8 Antibody

Catalog Number: A00042S347


Overview

Product Name	Anti-Caspase-8 (S347) CASP8 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Caspase-8 (S347) CASP8 Antibody catalog # A00042S347. Tested in WB,IHC,IP applications. This antibody reacts with Human,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14790


Technical Details

Immunogen	Synthetic peptide, corresponding to Human Caspase-8.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IP: 1:50-1:200

Anti-Caspase-8 (S347) CASP8 Antibody (A00042S347) Images



Western blot (WB) analysis of Caspase-8 pAb at 1:2000 dilution Lane1:PC12 whole cell lysate(30ug) Lane2:BV2 whole cell lysate(30ug) Lane3:A549 whole cell lysate(30ug) Lane4:MCF-7 whole cell lysate(20ug)



Immunohistochemistry (IHC) analyzes of Caspase-8 (S347) pAb in paraffin-embedded human lung carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caspase-8 (S347) CASP8 Antibody

For Research Use Only. Not for use in diagnostic procedures.