

## Anti-KMT6/EZH2 Antibody

Catalog Number: A00050-2


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-KMT6/EZH2 Antibody  |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-KMT6/EZH2 Antibody catalog # A00050-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human.              |
| Application          | IF, IHC, ICC, WB   |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q15910   |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | recombinant protein                                       |
| Cross Reactivity    | No cross reactivity with other proteins.                  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1.0 mg/mL by UV absorbance at 280 nm                      |
| Purification        | ProA affinity purified                                    |
| Suggested Dilutions | WB: 1:1,000-1:2,000<br>ICC: 1:50-1:200<br>IHC: 1:50-1:200 |

## Anti-KMT6/EZH2 Antibody (A00050-2) Images



Western blot analysis of KMT6/EZH2 on different lysates using anti-KMT6/EZH2 antibody at 1/1,000 dilution. Positive control: Lane 1: 293T Lane 2: 293



ICC staining KMT6/EZH2 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KMT6/EZH2 Antibody

For Research Use Only. Not for use in diagnostic procedures.