

Anti-Smad3 (phospho-S204) Antibody

Catalog Number: A00059S204

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Smad3 (phospho-S204) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Smad3 (phospho-S204) Antibody catalog # A00059S204. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84022

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Smad3 around the phosphorylation site of Serine 204.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 100% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Smad3 (phospho-S204) Antibody (A00059S204) Images



Anti-Smad3 (phospho-S204) Antibody

For Research Use Only. Not for use in diagnostic procedures.