

## Anti-IGF-1R (phospho-Y980) Antibody

Catalog Number: A00070Y980

### Overview


Product Name	Anti-IGF-1R (phospho-Y980) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IGF-1R (phospho-Y980) Antibody catalog # A00070Y980. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08069

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human IGF-1R around the phosphorylation site of Tyrosine 980.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-IGF-1R (phospho-Y980) Antibody (A00070Y980) Images

---



Immunohistochemistry (IHC) analyzes of p-IGF1R (Y980) pAb in paraffin-embedded human lung adenocarcinoma tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IGF-1R (phospho-Y980) Antibody

For Research Use Only. Not for use in diagnostic procedures.