

Anti-AR-Beta 2 ADRB2 Antibody

Catalog Number: A00072-1


Overview

Product Name	Anti-AR-Beta 2 ADRB2 Antibody
Reactive Species	Human, Mouse, Rat, Zebrafish
Description	Boster Bio Anti-AR-Beta 2 ADRB2 Antibody catalog # A00072-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat,zebrafish.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07550


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 IHC: 1:50-1:200

Anti-AR-Beta 2 ADRB2 Antibody (A00072-1) Images



Western blot analysis of ADRB2 on zebrafish lysates using anti-ADRB2 antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded mouse liver tissue using anti-ADRB2 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AR-Beta 2 ADRB2 Antibody

For Research Use Only. Not for use in diagnostic procedures.