

## Anti-Raf-B (A597) BRAF Antibody

Catalog Number: A00075-1

### Overview

Product Name	Anti-Raf-B (A597) BRAF Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Raf-B (A597) BRAF Antibody catalog # A00075-1. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15056

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 560-610 of Human Raf-B.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-Raf-B (A597) BRAF Antibody (A00075-1) Images

---




Figure 1. Immunohistochemistry validation of BRAF using Anti-Raf-B (A597) BRAF Antibody (A00075-1).

Immunohistochemistry (IHC) analyzes of Raf-B (A597) pAb in paraffin-embedded human breast carcinoma tissue at 1:100. For more protocol information of IHC

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Raf-B (A597) BRAF Antibody