

Anti-B-Raf (Phospho-T401) Antibody

Catalog Number: A00075T401


Overview

Product Name	Anti-B-Raf (Phospho-T401) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-B-Raf (Phospho-T401) Antibody catalog # A00075T401. Tested in WB, ICC, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15056


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr401 of human B Raf.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:100-1:500

Anti-B-Raf (Phospho-T401) Antibody (A00075T401) Images



Western blot analysis of Phospho-B Raf(T401) on different lysates using anti-Phospho-B Raf(T401) antibody at 1/1,000 dilution. Positive control: Lane 1: Mouse testis Lane 2: PC-12



ICC staining Phospho-B Raf(T401) in PC-3M cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-B-Raf (Phospho-T401) Antibody

For Research Use Only. Not for use in diagnostic procedures.