

Anti-Raf-B (Y198) BRAF Antibody

Catalog Number: A00075Y198


Overview

Product Name	Anti-Raf-B (Y198) BRAF Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Raf-B (Y198) BRAF Antibody catalog # A00075Y198. Tested in IHC,IF applications. This antibody reacts with Human,Mouse.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15056

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 170-220 of Human Raf-B.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Raf-B (Y198) BRAF Antibody (A00075Y198) Images



Immunohistochemistry (IHC) analyzes of Raf-B (Y198) pAb in paraffin-embedded human brain carcinoma tissue at 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Raf-B (Y198) BRAF Antibody

For Research Use Only. Not for use in diagnostic procedures.