

## Anti-Interleukin-13 IL13 Antibody

Catalog Number: A00077-1


### Overview

Product Name	Anti-Interleukin-13 IL13 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Interleukin-13 IL13 Antibody catalog # A00077-1. Tested in IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35225

### Technical Details

Immunogen	Recombinant full length Human IL13.
Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:2000 IHC: 1:50~1:200</p> <p>For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a></p>

### Anti-Interleukin-13 IL13 Antibody (A00077-1) Images



Western blot analysis of extracts of recombinant protein using IL13 antibody.

## 1 Publications Citing This Product

1. PubMed ID: 27916586, Selaginella uncinata flavonoids ameliorated ovalbumin-induced airway inflammation in a rat model of asthma

Visit [bosterbio.com/anti-interleukin-13-il13-antibody-a00077-1-boster.html](https://bosterbio.com/anti-interleukin-13-il13-antibody-a00077-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-13 IL13 Antibody