

Anti-Caspase 9 (Tyr153) CASP9 Antibody

Catalog Number: A00080Y153


Overview

Product Name	Anti-Caspase 9 (Tyr153) CASP9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Caspase 9 (Tyr153) CASP9 Antibody catalog # A00080Y153. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P55211


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 130-180 of Human Caspase 9.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Caspase 9 (Tyr153) CASP9 Antibody (A00080Y153) Images



Western blot (WB) analysis of Caspase 9 (Tyr153) polyclonal antibody at 1:500 dilution Lane1:LO2 whole cell lysate(40ug) Lane2:LO2 treated with PBS(1×PBS,PH7.4) for 30 minutes whole cell lysate(40ug) Lane3:LO2 treated with PBS(1×PBS,PH7.4) for 120 minutes whole cell lysate(40ug) Lane4:LO2 treated with PBS(1×PBS,PH7.4) for 240 minutes whole cell lysate(40ug)



Immunofluorescence analysis of Jurkat cells using Caspase 9 (Tyr153) pAb at dilution of 1:50 (40x lens).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caspase 9 (Tyr153) CASP9 Antibody

For Research Use Only. Not for use in diagnostic procedures.