

Anti-Rad51 (T309) Antibody

Catalog Number: A00088T309


Overview

Product Name	Anti-Rad51 (T309) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Rad51 (T309) Antibody catalog # A00088T309. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q06609

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 280-330 of Human Rad51.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Rad51 (T309) Antibody (A00088T309) Images



Immunohistochemistry (IHC) analyzes of Rad51 (T309) pAb in paraffin-embedded human colon carcinoma tissue.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rad51 (T309) Antibody

For Research Use Only. Not for use in diagnostic procedures.