

Anti-Smad2/3 (Phospho-T8) Antibody

Catalog Number: A00090T8


Overview

Product Name	Anti-Smad2/3 (Phospho-T8) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Smad2/3 (Phospho-T8) Antibody catalog # A00090T8. Tested in WB applications. This antibody reacts with Human,Rat,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15796/P84022

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Smad2/3 around the phosphorylation site of Threonine 8.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Smad2/3 (Phospho-T8) Antibody (A00090T8) Images



Western blot (WB) analysis of Smad2/3 (Phospho-T8) polyclonal antibody at 1:500 dilution Lane1:C6 whole cell lysate(40ug) Lane2:MEF whole cell lysate(40ug) Lane3:A2780 whole cell lysate(40ug) Lane4:A549 whole cell lysate(40ug) Lane5:A375 whole cell lysate(40ug) Lane6:MCF-7 whole cell lysate(30ug) Lane7:Hela whole cell lysate(30ug)

1 Publications Citing This Product

1. PubMed ID: 34093874, Wang Q,Liu J,Hu Y,Pan T,Xu Y,Yu J,Xiong W,Zhou Q,Wang Y.Local administration of liposomal-based SrpX2 gene therapy reverses pulmonary fibrosis by blockading fibroblast-to-myofibroblast transition.Theranostics.2021 May 13;11(14):7110-7125.doi:10.7150/thno.61085. PMID:34093874;PMCID:PMC8171094.

Visit bosterbio.com/anti-smad2-3-phospho-t8-antibody-a00090t8-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Smad2/3 (Phospho-T8) Antibody

For Research Use Only. Not for use in diagnostic procedures.