

Anti-PDGFR Beta (phospho Y740) Pdgfrb Antibody

Catalog Number: A00096Y740


Overview

Product Name	Anti-PDGFR Beta (phospho Y740) Pdgfrb Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-PDGFR Beta (phospho Y740) Pdgfrb Antibody catalog # A00096Y740. Tested in WB, ICC/IF applications. This antibody reacts with Mouse, Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05622


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Tyr740 of human PDGF Receptor beta.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC/IF: 1:50-1:200

Anti-PDGFR Beta (phospho Y740) Pdgfrb Antibody (A00096Y740) Images



ICC staining PDGF Receptor beta(phospho Y740) in C2C12 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining PDGF Receptor beta(phospho Y740) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDGFR Beta (phospho Y740) Pdgfrb Antibody

For Research Use Only. Not for use in diagnostic procedures.