

Anti-HLA-DQB1/B2 Antibody

Catalog Number: A00106-2

Overview


Product Name	Anti-HLA-DQB1/B2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HLA-DQB1/B2 Antibody catalog # A00106-2. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01920

Technical Details

Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:1000 For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/

Anti-HLA-DQB1/B2 Antibody (A00106-2) Images

Figure 1. Western blotting validation for Anti-HLA-DQB1/B2 Antibody A00106-2



Western blot (WB) analysis of HLA-DQB1/B2 polyclonal antibody at 1:500 dilution
Lane1:HEK293T whole cell lysate
Lane2:Raw264.7 whole cell lysate
Lane3:PC12 whole cell lysate
Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HLA-DQB1/B2 Antibody