

## Anti-Ribonuclease 3 Drosha Antibody

Catalog Number: A00111-2


### Overview

Product Name	Anti-Ribonuclease 3 Drosha Antibody
Reactive Species	Human
Description	Boster Bio Anti-Ribonuclease 3 Drosha Antibody catalog # A00111-2. Tested in WB,ICC,IF,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NRR4


### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

## Anti-Ribonuclease 3 Drosha Antibody (A00111-2) Images



Western blot analysis of Drosha on HeLa cell and SiHa cell lysates using anti-Drosha antibody at 1/500 dilution.



ICC staining Drosha in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ribonuclease 3 Drosha Antibody

For Research Use Only. Not for use in diagnostic procedures.