

Anti-Ras (H27) HRAS Antibody

Catalog Number: A00114-1


Overview

Product Name	Anti-Ras (H27) HRAS Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ras (H27) HRAS Antibody catalog # A00114-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01112/P01116/P01111


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 1-50 of Human Ras.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Ras (H27) HRAS Antibody (A00114-1) Images



Western blot (WB) analysis of Ras (H27) polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate Lane2:Hela cell lysate Lane3:Mouse brain tissue lysate Lane4:Rat heart tissue lysate



Western blot (WB) analysis of Ras (H27) pAb at 1:500 dilution Lane1:Hela whole cell lysate(40ug) Lane2:NIH-3T3 whole cell lysate(40ug) Lane3:PC12 whole cell lysate(40ug)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ras (H27) HRAS Antibody

For Research Use Only. Not for use in diagnostic procedures.