

Anti-YAP1 (Phospho-S127) Antibody

Catalog Number: A00116S127


Overview

Product Name	Anti-YAP1 (Phospho-S127) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-YAP1 (Phospho-S127) Antibody catalog # A00116S127. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P46937


Technical Details

Immunogen	Human breast carcinoma tissue, human kidney tissue, mouse kidney tissue.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

Anti-YAP1 (Phospho-S127) Antibody (A00116S127) Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-phospho-YAP1 (S127) antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-phospho-YAP1 (S127) antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-YAP1 (Phospho-S127) Antibody

For Research Use Only. Not for use in diagnostic procedures.