

Anti-IKK Beta (phospho-Y199) IKBKB Antibody

Catalog Number: A00118Y199


Overview

Product Name	Anti-IKK Beta (phospho-Y199) IKBKB Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IKK Beta (phospho-Y199) IKBKB Antibody catalog # A00118Y199. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14920

Technical Details

Immunogen	Synthetic phosphopeptide derived from human eEF2K around the phosphorylation site of Tyrosine 199.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-IKK Beta (phospho-Y199) IKBKB Antibody (A00118Y199) Images



Immunohistochemistry (IHC) analyzes of p-IKKbeta (Y199) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IKK Beta (phospho-Y199) IKBKB Antibody

For Research Use Only. Not for use in diagnostic procedures.