

Anti-Cystathionine beta-synthase CBS Antibody

Catalog Number: A00130


Overview

Product Name	Anti-Cystathionine beta-synthase CBS Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cystathionine beta-synthase CBS Antibody catalog # A00130. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35520


Technical Details

Immunogen	Recombinant fusion protein of human CBS(NP_000062.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200

Anti-Cystathionine beta-synthase CBS Antibody (A00130) Images



Western blot analysis of CBS expression in HeLa (A), MCF7 (B), Romas (C), mouse brain (D) whole cell lysates.



Immunofluorescent analysis of CBS staining in U2OS cells.

1 Publications Citing This Product

1. PubMed ID: 26798657, Heme Oxygenase-1 Promotes Delayed Wound Healing in Diabetic Rats

Visit bosterbio.com/anti-cystathionine-beta-synthase-cbs-antibody-a00130-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cystathionine beta-synthase CBS Antibody

For Research Use Only. Not for use in diagnostic procedures.