

Anti-Huntingtin HTT Antibody

Catalog Number: A00134-1

Overview

Product Name	Anti-Huntingtin HTT Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Huntingtin HTT Antibody catalog # A00134-1. Tested in WB,IHC,Flow Cytometry,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42858


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50-1:200 FC: 1:50-1:100 ICC: 1:50

Anti-Huntingtin HTT Antibody (A00134-1) Images



Western blot analysis of Huntingtin on A431 cell using anti-Huntingtin antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue using anti-Huntingtin antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Huntingtin HTT Antibody

For Research Use Only. Not for use in diagnostic procedures.