

Anti-DNMT3A Antibody

Catalog Number: A00136


Overview

Product Name	Anti-DNMT3A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DNMT3A Antibody catalog # A00136. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y6K1

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human DNMT3A
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-DNMT3A Antibody (A00136) Images



Western blot (WB) analysis of DNMT3A pAb at 1:500 dilution
Lane1:MEF whole cell lysate(40ug) Lane2:H9C2 whole cell lysate(40ug) Lane3:HCT116 whole cell lysate(40ug)
Lane4:Jurkat whole cell lysate(40ug) Lane5:PC3 whole cell lysate (40ug) Lane6:HEK293T whole cell lysate (40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNMT3A Antibody

For Research Use Only. Not for use in diagnostic procedures.