

Anti-p21 (phospho-T145) CDKN1A Antibody

Catalog Number: A00145T145


Overview

Product Name	Anti-p21 (phospho-T145) CDKN1A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p21 (phospho-T145) CDKN1A Antibody catalog # A00145T145. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P38936

Technical Details

Immunogen	Synthetic phosphopeptide derived from human p21 around the phosphorylation site of Threonine 145.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-p21 (phospho-T145) CDKN1A Antibody (A00145T145) Images



Immunohistochemistry (IHC) analyzes of p-p21 (T145) pAb in paraffin-embedded human tonsil carcinoma tissue at 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p21 (phospho-T145) CDKN1A Antibody

For Research Use Only. Not for use in diagnostic procedures.