

## Anti-FAK (Y570) PTK2 Antibody

Catalog Number: A00151Y570

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-FAK (Y570) PTK2 Antibody  |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-FAK (Y570) PTK2 Antibody catalog # A00151Y570. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.   |
| Application          | IHC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q05397   |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthetic peptide, corresponding to amino acids 550-600 of Human FAK.   |
| Cross Reactivity    | No cross reactivity with other proteins.  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000<br>IHC: 1:50-1:200   |

## Anti-FAK (Y570) PTK2 Antibody (A00151Y570) Images



Western blot (WB) analysis of FAK (Y570) polyclonal antibody at 1:500 dilution Lane1:HCT116 whole cell lysate(40ug) Lane2:LO2 whole cell lysate(40ug) Lane3:HepG2 whole cell lysate(40ug) Lane4:BV2 whole cell lysate(40ug) Lane5:C6 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of FAK (Y570) pAb in paraffin-embedded human brain tissue.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FAK (Y570) PTK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.