

## Anti-Cyclin-dependent kinase 4 Cdk4 Antibody

Catalog Number: A00159-1


### Overview

Product Name	Anti-Cyclin-dependent kinase 4 Cdk4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cyclin-dependent kinase 4 Cdk4 Antibody catalog # A00159-1. Tested in WB,ICC/IF,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11802


### Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:100-1:500 IHC: 1:50-1:200 FC: 1:50-1:100

## Anti-Cyclin-dependent kinase 4 Cdk4 Antibody (A00159-1) Images



Western blot analysis of Cdk4 on HeLa cells lysates using anti-Cdk4 antibody at 1/1,000 dilution.



ICC staining Cdk4 in HeLa cells (green). The nuclear counterstain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin-dependent kinase 4 Cdk4 Antibody

For Research Use Only. Not for use in diagnostic procedures.