

Anti-Interferon regulatory factor 3 IRF3 Antibody

Catalog Number: A00165-1

About IRF3

Interferons (IFN)s are involved in a multitude of immune interactions during viral infections and play a major role in both the induction and regulation of innate and adaptive antiviral mechanisms. During infection, host-virus interactions signal downstream molecules such as transcription factors such as IFN regulatory factor-3 (IRF3) which can act to stimulate transcription of IFN-alpha/beta genes. IRF3 is present in an inactive form in the cytoplasm of most cells. Following viral infection, IRF3 can be activated by I κ B kinase-epsilon and TANK-binding kinase 1 (TBK1), whereupon IRF3 translocates to the nucleus. IRF3 can also be activated by stimulation of toll-like receptor 3 (TLR3) by dsRNA. IRF3 exists as at least two distinct isoforms.

Overview

Product Name	Anti-Interferon regulatory factor 3 IRF3 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Interferon regulatory factor 3 IRF3 Antibody (Catalog # A00165-1). Tested in ELISA, WB, ICC, IF applications. This antibody reacts with Human, Rat.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	IRF3 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	IRF3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q14653


Technical Details

Immunogen	IRF3 antibody antibody was raised against a peptide corresponding to 14 amino acids near the center of human IRF3. The immunogen is located within amino acids 150 - 200 of IRF3.
Predicted Reactive Species	Bovine, Pig
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	IRF3 Antibody is affinity chromatography purified via peptide column.


Suggested Dilutions

IRF3 antibody can be used for detection of IRF3 by Western blot at 1 ug/mL. Antibody can also be used for immunocytochemistry starting at 2 ug/mL. For immunofluorescence start at 2 ug/mL. Antibody validated: Western Blot in human and rat samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-Interferon regulatory factor 3 IRF3 Antibody (A00165-1) Images



Western blot analysis of IRF3 in (A) human kidney and (B) rat kidney lysate with IRF3 antibody at 1 ug/mL.



Immunocytochemistry of IRF3 in Ramos cells with IRF3 antibody at 2 ug/mL.



Immunofluorescence of IRF3 in Ramos cells with IRF3 antibody at 2 ug/mL.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interferon regulatory factor 3 IRF3 Antibody

For Research Use Only. Not for use in diagnostic procedures.