

Anti-IRF3 (Phospho-S386) Antibody

Catalog Number: A00165S386


Overview

Product Name	Anti-IRF3 (Phospho-S386) Antibody
Reactive Species	Human
Description	Boster Bio Anti-IRF3 (Phospho-S386) Antibody catalog # A00165S386. Tested in WB,ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14653

Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser386 of human IRF3.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:50-1:200

Anti-IRF3 (Phospho-S386) Antibody (A00165S386) Images



Western blot analysis of IRF3 (Phospho-S386) expression in HeLa-LPS (A), HepG2 (B) whole cell lysates.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IRF3 (Phospho-S386) Antibody

For Research Use Only. Not for use in diagnostic procedures.