

Anti-Dnmt1 Antibody

Catalog Number: A00172


Overview

Product Name	Anti-Dnmt1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Dnmt1 Antibody catalog # A00172. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P26358

Technical Details

Immunogen	Recombinant fusion protein of human DNMT1(NP_001370.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200 IP: 1:50-1:100

Anti-Dnmt1 Antibody (A00172) Images



Immunoprecipitation analysis of 200ug extracts of 293T cells using 3ug DNMT1 antibody . Western blot was performed from the immunoprecipitate using DNMT1 antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dnmt1 Antibody

For Research Use Only. Not for use in diagnostic procedures.