

Anti-Dnmt1 Antibody

Catalog Number: A00172-2

Overview

Product Name	Anti-Dnmt1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Dnmt1 Antibody catalog # A00172-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P26358

Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-Dnmt1 Antibody (A00172-2) Images




Figure 2. Immunocytochemistry staining of DNMT1 using Anti-Dnmt1 Antibody (A00172-2).

ICC staining Dnmt1 in HeLa cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 1. Western blotting validation for Anti-Dnmt1 Antibody A00172-2

Western blot analysis of Dnmt1 on HepG2 cells lysates using anti-Dnmt1 antibody at 1/1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dnmt1 Antibody