

Anti-p27 Kip1 (phospho-T198) CDKN1B Antibody

Catalog Number: A00173T198


Overview

Product Name	Anti-p27 Kip1 (phospho-T198) CDKN1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p27 Kip1 (phospho-T198) CDKN1B Antibody catalog # A00173T198. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P46527

Technical Details

Immunogen	Synthetic phosphopeptide derived from human p27 Kip1 around the phosphorylation site of Threonine 198.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-p27 Kip1 (phospho-T198) CDKN1B Antibody (A00173T198) Images



Western blot (WB) analysis of p-p27 Kip1 (T198) polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate treated with EGF(0.1ng/ML,30mins) Lane2:Raw264.7 cell lysate treated with EGF(0.1ng/ML,30mins) Lane3:H9C2 cell lysate treated with EGF(0.1ng/ML,30mins)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p27 Kip1 (phospho-T198) CDKN1B Antibody

For Research Use Only. Not for use in diagnostic procedures.